




SOUTHERNIMPLANTS®
Innovative Treatment Solutions

SI-BASE Abutments



Abutment versatility redefined




SI-BASE ABUTMENTS

THE NEW GENERATION OF TITANIUM ABUTMENTS

Southern's SI-BASE range adds versatility to the clinician's armamentarium by providing a full range of customisable titanium abutments.

What variations of SI-BASE are available?




8 mm tall chimney

Laser-marked at 4.5 and 6 mm to enable easy modification of the chimney to match the crown height.



Yellow anodised for aesthetics

Providing less abutment shine through in thinner restorations.



Engaging and non-engaging

Single and multi-unit restorations.

Collar heights

3 different collar height variations enabling ideal tissue profiling for better emergence.
1.5 mm, 3 mm and 5 mm (additional 0 mm in the IT range only).

Platform diameters

2 restorative platform widths:
Narrow (in collar heights 0 mm, 1.5 mm and 3 mm) and wide (in collar heights 3 mm and 5 mm) variations to better support crown emergence.


SI-Digital workflow

The SI-BASE abutments are compatible with digital workflows in both



Efficiency

Due to multiple variations, clinicians can select the ideal abutment requiring the least amount of modification.



SI-BASE Abutments

External Hex

Implant connection interface diameter	Collar height (mm)	Restorative platform		Retaining screw	Indications
		Narrow	Wide		
Ø3.0 (Piccolo (IP))	1.5	SIB-EX-30-40C1.5 SIB-NX-30-40C1.5		TS-P-16	Maxillary and mandibular lateral incisors and mandibular central incisor
	3	SIB-EX-30-40C3 SIB-NX-30-40C3	SIB-EX-30-45C3 SIB-NX-30-45C3		
	5		SIB-EX-30-45C5 SIB-NX-30-45C5		
Ø3.25 (IBN)	1.5	SIB-EX-34-43C1.5 SIB-NX-34-43C1.5		3 Series screws	Maxillary and mandibular lateral and central incisors
	3	SIB-EX-34-43C3 SIB-NX-34-43C3	SIB-EX-34-48C3 SIB-NX-34-48C3		
	5		SIB-EX-34-48C5 SIB-NX-34-48C5		
Ø4.0 (IB)	1.5	SIB-EX-40-50C1.5 SIB-NX-40-50C1.5		3 Series screws	All regions of the mouth
	3	SIB-EX-40-50C3 SIB-NX-40-50C3	SIB-EX-40-55C3 SIB-NX-40-55C3		
	5		SIB-EX-40-55C5 SIB-NX-40-55C5		
Ø5.0 (BA)	1.5	SIB-EX-50-55C1.5 SIB-NX-50-55C1.5		3 Series screws	All regions of the mouth
	3	SIB-EX-50-55C3 SIB-NX-50-55C3	SIB-EX-50-65C3 SIB-NX-50-65C3		
	5		SIB-EX-50-65C5 SIB-NX-50-65C5		
Ø6.0 (BBB)	1.5	SIB-EX-60-65C1.5 SIB-NX-60-65C1.5		3 Series screws	All regions of the mouth
	3	SIB-EX-60-65C3 SIB-NX-60-65C3	SIB-EX-60-75C3 SIB-NX-60-75C3		
	5		SIB-EX-60-75C5 SIB-NX-60-75C5		
Ø7.0 (MAX 9)	1.5	SIB-EX-70-75C1.5 SIB-NX-70-75C1.5		3 Series screws	All regions of the mouth
	3	SIB-EX-70-75C3 SIB-NX-70-75C3	SIB-EX-70-80C3 SIB-NX-70-80C3		
	5		SIB-EX-70-80C5 SIB-NX-70-80C5		

Internal Octagon (IT)

Implant connection interface diameter	Collar height (mm)	Restorative platform		Retaining screw	Indications
		Narrow	Wide		
ITC3 ITC4/ITT4 ITC5/ITT5 MAXIT7	0	SIB-ITS-48 SIB-ITSNE-48		TSIT2	All regions of the mouth
	1.5	SIB-ITS-55C1.5 SIB-ITSNE-55C1.5	SIB-ITS-60C1.5 SIB-ITSNE-60C1.5		
	3	SIB-ITS-55C3 SIB-ITSNE-55C3			
ITC6-5/ITT6-5 ITT6 MAXIT8 MAXIT9	0	SIB-IT6-65 SIB-IT6NE-65		TSIT2	All regions of the mouth
	1.5		SIB-IT6-75C1.5 SIB-IT6NE-75C1.5		

SI-BASE Abutments

TRI-NEX®

Implant connection interface diameter	Collar height (mm)	Restorative platform		Retaining screw	Indications
		Narrow	Wide		
TRI-NEX 35	1.5	SIB-EL-35-43C1.5 SIB-NL-35-43C1.5		TS-L-18	All regions of the mouth
	3	SIB-EL-35-43C3 SIB-NL-35-43C3	SIB-EL-35-48C3 SIB-NL-35-48C3		
	5		SIB-EL-35-48C5 SIB-NL-35-48C5		
TRI-NEX 43	1.5	SIB-EL-43-45C1.5 SIB-NL-43-45C1.5		TS-L-20	All regions of the mouth
	3	SIB-EL-43-45C3 SIB-NL-43-45C3	SIB-EL-43-55C3 SIB-NL-43-55C3		
	5		SIB-EL-43-55C5 SIB-NL-43-55C5		
TRI-NEX 50	1.5	SIB-EL-50-55C1.5 SIB-NL-50-55C1.5		TS-L-20	All regions of the mouth
	3	SIB-EL-50-55C3 SIB-NL-50-55C3	SIB-EL-50-65C3 SIB-NL-50-65C3		
	5		SIB-EL-50-65C5 SIB-NL-50-65C5		
TRI-NEX 60	1.5	SIB-EL-60-65C1.5 SIB-NL-60-65C1.5		TS-L-20	All regions of the mouth
	3	SIB-EL-60-65C3 SIB-NL-60-65C3	SIB-EL-60-75C3 SIB-NL-60-75C3		

Deep Conical (DC)

Implant connection interface diameter	Collar height (mm)	Restorative platform		Retaining screw	Indications
		Narrow	Wide		
DC3	1.5	SIB-DC3-40C1.5 SIB-NDC3-40C1.5		TS-DC3-14	Maxillary and mandibular lateral incisors and mandibular central incisor
	3	SIB-DC3-40C3 SIB-NDC3-40C3	SIB-DC3-45C3 SIB-NDC3-45C3		
	5		SIB-DC3-45C5 SIB-NDC3-45C5		
DC4	1.5	SIB-DC4-43C1.5 SIB-NDC4-43C1.5		TS-DC4-16	All regions of the mouth
	3	SIB-DC4-43C3 SIB-NDC4-43C3	SIB-DC4-48C3 SIB-NDC4-48C3		
	5		SIB-DC4-48C5 SIB-NDC4-48C5		
DC5	1.5	SIB-DC5-55C1.5 SIB-NDC5-55C1.5		TS-DC5-20	All regions of the mouth
	3	SIB-DC5-55C3 SIB-NDC5-55C3	SIB-DC5-65C3 SIB-NDC5-65C3		
	5		SIB-DC5-65C5 SIB-NDC5-65C5		

Compact Conical Abutments (TMC)

Abutment connection interface diameter	Collar height (mm)	Restorative platform		Retaining screw	Indications
		Item code	Platform diameter		
Ø4.8	1	SIB-TMC1	Ø5.2	1 Series screw	All regions of the mouth
Ø6.0	1	SIB-TMCW1	Ø6.4	1 Series screw	All regions of the mouth

SI-BASE Abutments

Internal Hex (M-Series & PROVATA®)

Implant connection interface diameter	Collar height (mm)	Restorative platform		Retaining screw	Indications
		Narrow	Wide		
PRO3	1.5	SIB-3M-40C1.5 SIB-3NM-40C1.5		TS-Z-16	Maxillary and mandibular lateral and central incisors
	3	SIB-3M-40C3 SIB-3NM-40C3	SIB-3M-45C3 SIB-3NM-45C3		
	5		SIB-3M-45C5 SIB-3NM-45C5		
PRO4 PRO5 PROMAX6 IM-T37 IM-T43 IM-T50	1.5	SIB-M-43C1.5 SIB-NM-43C1.5		TS-Z-18	All regions of the mouth
	3	SIB-M-43C3 SIB-NM-43C3	SIB-M-48C3 SIB-NM-48C3		
	5		SIB-M-48C5 SIB-NM-48C5		
PRO5 PROMAX6 IM-T50 (platform matched)	1.5	SIB-M-PM-48C1.5		TS-Z-18	All regions of the mouth
	3	SIB-M-PM-48C3	SIB-M-PM-55C3		
	5		SIB-M-PM-55C5		
MSC-PRO6 PROMAX7/8/9	1.5	SIB-Z-65C1.5 SIB-NZ-65C1.5		TS-Z-18	All regions of the mouth
	3	SIB-Z-65C3 SIB-NZ-65C3	SIB-Z-75C3 SIB-NZ-75C3		
	5		SIB-Z-75C5 SIB-NZ-75C5		
PROMAX8/9 (platform matched)	3		SIB-Z-PM-75C3	TS-Z-18	All regions of the mouth
	5		SIB-Z-PM-75C5		

Single Platform

Implant connection interface diameter	Collar height (mm)	Restorative platform		Retaining screw	Indications
		Narrow	Wide		
SP35 SP40 SP50	1.5	SIB-SP-43C1.5 SIB-NSP-43C1.5		TS-SP-16	All regions of the mouth
	3	SIB-SP-43C3 SIB-NSP-43C3	SIB-SP-48C3 SIB-NSP-48C3		
	5		SIB-SP-48C5 SIB-NSP-48C5		
SP50 (platform matched)	1.5	SIB-SP-55-PM-C1.5		TS-SP-16	All regions of the mouth
	3	SIB-SP-55-PM-C3	SIB-SP-65-PM-C3		
	5		SIB-SP-65-PM-C5		



South Africa - Headquarters 1 Albert Road, Irene, RSA
T: +27-12-667-1046 | E: info@southernimplants.com

EC REP

Southern Implants Europe AB Holmgatan 30, S-791 71, Falun, Sweden
E: ecrep@southernimplants.com

Subsidiaries

Australia
Southern Implants Australia
T: +61-2-8076-9337
E: info@southernimplants.com.au

Spain and Portugal
Southern Implants Ibérica
T: +34 935 053 507
E: info@southernimplants.com.es

United Kingdom and Ireland
Southern Implants UK
T: +44-20-8899-6845 / 6 / 7
E: info@southernimplants.com.uk

USA and Canada
Southern Implants North America Inc.
T: +1-561-472-0990
E: customercare@southernimplants.com

Southern Implants® are distributed world-wide, please visit southernimplants.com for a list of Distributors.
For the latest and full assortment of Southern Implant products, visit our current resources library at southernimplants.com.

CAT-6008-00